





Large screen smart touch control



Small-caliber nozzle



Automatic identification of scaling mode



Multi-function foot pedal



Elliptical vibration

The elliptical vibration trajectory can achieve 0 ° angle operation between the tip and the tooth surface, enabling the polishing effect.

Disperse the impact force between the tip and the teeth to achieve comfortable scaling. Achieve 360° scaling.

Titanium alloy tip

The hardness is lower than dental enamel, more flexible and more comfortable.

Titanium alloy tip hardness: HRC30 Stainless steel hardness: HRC50



Implant maintenance





Sales department: 0773-5873196/2350599 Fax: 0773-5822450

E-mail: woodpecker4@glwoodpecker.com Website: http://www.glzmn.com, http://www.glwoodpecker


texDAN

texDAN s.r.o., Chelčického 1143, 542 32 Úpice; IČ: 287 68 591 mobil: +420 739 693 258, +420 603 558 500 info@woodpecker.cz www.woodpecker.cz

Air polishing is the process of using compressed air to spray powder and water on the tooth surface with Air Polisher to remove tooth surface adhesive substances (such as dental plaque, pigment, soft dirt, etc.) and polish the tooth surface. Efficient and comfortable treatment makes it more endurable for patients. Does not scratch the surface of the teeth and does not easily damage soft tissues.

Effective protection of teeth surface

The hardness of the powder is much lower than that of the tooth enamel, so the scaling treatment can better protect the teeth, avoiding tooth surface damage caused by traditional over-scraping, and making the treatment more endurable for patients.



Thorough removal of plaque

Air polishing can achieve significant scaling effect in the occlusal surface, adjacent surface, orthodontic brackets and other difficult-to-clean areas.

Reduce plaque and pigment adhesion

After ultrasonic scaling or manual scaling, patients are likely to feel that the teeth surface is rough and uneven. The main function of air polishing is to polish teeth surface. After polishing, the tooth surface is smoother, and the plaque and pigment are less likely to adhere to it.



Plaque Control

Small-caliber nozzle, more concentrated air polishing, and more efficient plaque removal High concentricity of water and air path achieves a good coating of powder by water and reduces the



Scaling effect



Dual water supply

Equipped with two water supply bottles (600mL, 1400mL) made of anti-corrosive materials. They can contain special chemicals such as chlorhexidine, hydrogen peroxide, and sodium hypochlorite. External water supply is also available.

Water heating unit makes the treatment more endurable.

Three-section Air polishing Handpiece, easier for cleaning

Lightweight handpiece Longer service life



Classic sodium bicarbonate powder

Average particle size: about 65 microns; hardness is less than that of enamel, doing no harm to teeth Contain natural plant-derived essence for a more pleasant


The powder is waterproof to prevent the device from being blocked by wet powder.

Teeth	Moh's hardness	Powder	Moh's hardne
Enamel	4-5	Sodium bicarbonate	2.5
Dentin	2-2.5	Glycine	2
		Calcium carbonate	3
Cementum	2-4	Aluminum hydroxide	3
•	•		•

Funnel-shaped powder tank, more stable powder flow

It can be quickly plugged and removed for easily Working Time and Powder Consumption of PT-A





Multi-function foot pedal








Root canal irrigation







with security, accuracy and

Ultrasonic periodontal treatment


Save supragingival and subgingival ultrasonic scaling.

Elliptical vibration track disperses the impact force between the tip and teeth. The hardness of titanium alloy tip is lower than that of cementum and enamel, doing no harm to teeth.






Under scaling



Pick up the Ultrasonic Handpiece to automatically enter the Ultrasonic Scaling Mode



Pick up the Air Polishing Handpiece to automatically enter the Air Polishing Mode

Accessories

Main unit + Air-powder sterilizing box and Ultrasonic sterilizing box

Air-powder sterilizing box Polishing Handpiece and cleaning needle

For storage of all PT-A accessories

Standard configuration

Tip P59

Tip P33

Ultrasonic sterilizing boxFor sterilization and storage of Ultrasonic Handpiece and cleaning needle

Ultrasonic handpiece ×1 Woodpecker powder ×2 Titanium alloy tip P56×2、P59×1、P50L×2、P50R×2 Implant maintenance tip P94×1、P95×1

Endo irrigation tip E60×2

×1 | Power adapter | ×1