

**Valid from January 15th to June 28th 2024**



Only for participating dealers.  
All prices are in €, exclusive VAT.

**HUB Country Manager****Marc COLIN**

Tel.: +32 (0) 496 529 763

marc.colin@acteongroup.com

**Area Sales Manager Belux  
Equipment & Pharma****Raphaël NEF**

Tel.: +32 (0) 472 711 580

raphael.nef@acteongroup.com

**Area Sales Manager Imaging  
The Netherlands****Andy LUBERTI**

Tel.: +31 613 948 984

andy.luberti@acteongroup.com

**Area Sales Manager Surgery &  
Pharma The Netherlands****Sarah VAN ZANTEN**

Tel.: +31 630 815 978

sarah.vanzanten@acteongroup.com

**Area Sales Manager Imaging  
Nordic****Pär MÅRTENSSON**

Tel.: +46 (0) 738 473 733

par.martensson@acteongroup.com

**Area Sales Manager Nordic  
Equipment & Pharma****Marie ÅHGREN**

Tel.: +46 (0) 707 722 755

marie.ahgren@acteongroup.com

**Administration Commerciale****Sabine PARIS**

Belgium / Nordics

Tel.: +33 (0) 556 349 332

Fax: +33 (0) 566 349 331

sabin.paris@acteongroup.com

**Nirina INISAN**

Luxembourg / The Netherlands

Tel.: +33 (0) 556 349 250

Fax: +33 (0) 566 349 331

nirina.inisan@acteongroup.com

---

**ACTEON CUSTOMER SERVICE**

ACTEON's customer service  
can be contacted at:

**Imaging****PSPIX - SOPIX - CAMERAS - X MIND**

Tel.: +33 (0) 186 656 928

cs@acteongroup.com

**X MIND TRIUM - X MIND PRIME**

Tel.: +39 (0) 331 37 67 62

imaging.italysupport@acteongroup.com

**Equipment / Pharma**

Tel.: +33 (0) 556 349 163

cs@acteongroup.com



# Our Values

## **A human focused, service-oriented company**

The customer knows best: ACTEON®'s solutions are inspired and guided by the specific need of its diverse partners, be they healthcare professionals or distributors.

ACTEON®'s long-standing presence in the market has enabled the company to develop granular understanding of each type of customer and to create best-fit solutions.

## **Innovator at heart**

Creative thinking has always been a part of ACTEON®'s DNA and is still a core pillar of the company today. ACTEON® is building on this strength everyday to set new standards for the industry and improve people's lives.

## **Excellence and expertise**

ACTEON® has always stood for high quality and high performance. From advanced technologies to the manufacturing of manual instrumentation, the company has made a name for itself and grown trust through its requirements and its extensive experience in the industry.

## **Entrepreneurial spirit**

ACTEON® has created a work environment that enables collaborators to make the most of their capabilities, empowers them and encourages creativity and action.

## **Cooperation & knowledge spirit**

Team spirit is playing a central role in connecting the different entities of the group. It also infuses ACTEON®'s approach to its customers and partners and is a guiding principle in its ambition to transform the industry and everyday practice.

# Exclusive and patented Color Coding System™

Power selection is accurately determined with the Color Coding System™, which intuitively associates each tip color ring with one of the 4 power ranges







# newtron p5XS

~~€3.768,21~~  
**€2.099**

**P5XS b-LED** incl. with the following tips:  
Ref. F62100



# Newtron



Linear and regular vibrations

**Painless treatment for the patient**



Thin and light handpiece

**Decreased fatigue for the clinician**




No heat or handpiece vibration

**Better tactile sense**



Perfectly controlled vibrations

**Precise & gentle treatments**




Specific steel quality

**Anatomy & enamel preservation**




Total irrigation control

**No heat is safer for tissues**




Automatic and continuous frequency adjustments

**Maximal and constant efficiency whatever the load applied**



Real-time power compensation

**Continuous performances  
No force needed**



Powerful cavitation

**Deposits fragmentation  
Disinfecting effect**





**Newtron® handpiece**  
ref. F12281 - Autoclavable



~~€600,77~~  
**€315**

~~€831,57~~  
**€455**

**Newtron® LED handpiece**  
ref. F12609



**Newtron® Slim B. LED handpiece**  
> with blue LED ring: ref. F12900  
> with white LED ring: ref. F12905



~~€831,57~~  
**€509**

scaling 



**N°1 F00246**  
x1 €88 pce  
3+1 €264 pcs



**N°1S F00245**  
x1 €95 pce  
3+1 €285 pcs



**N°10P F00253**  
x1 €98 pce  
3+1 €294 pcs





**N°10Z F00254**  
x1 €129 pce  
3+1 €387 pcs



**N°10X F00359**  
x1 €98 pce  
3+1 €294 pcs

- > Excellent view, full LED ring
- > Avoid necrosis up to 120ml irrigation with peristaltic pump
- > Improve durability of the handpiece, external irrigation
- > Maximum hygiene, sterile single use or autoclavable lines
- > Work don't play, use of big bottles (1l) possible
- > Large choice of clinical application, 38 tips





**The alliance  
of technologies  
for powerful and secure  
bone surgeries**

### 3 modes:

- > Piezotome® (bone surgery).
- > Newtron® (perio, endo,...).
- > I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.

~~€14.183,77~~  
**€10.750**



### < INCLUSIVE



Ref. F12815 - Piezotome® LED handpiece + cord  
with a value of **€1592,18**  
Class IIa Medical devices (GMED) - CE0459



### Bone Surgery™ II Kit\*

Ref. F87509 - Tips BS1S II, BS2L II, BS2R II, BS4 II, BS5 II, BS6 II.

~~€1.356,85~~ **€1.090**



### Sinus Lift II Kit\*

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

~~€1.049,60~~ **€990**



### Intralift™ II Kit\*

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

~~€1.331,76~~ **€1.090**



### Extraction™ II Kit\*

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

~~€1.030,24~~ **€990**

\* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

# piezotome cube



~~€7.423,42~~  
**€5.650**



## Dynamic power adaptation for safe bone surgeries

Ref. F50100 - Delivered with a multifunction footswitch, a Piezotome® CUBE LED handpiece, 1 handpiece holder, an Essential II kit (6 tips BS1S, BS4, SL1, SL2, SL3 et LC2), a Piezotome® dynamometric wrench, 1 Irrigation pole, 2 autoclavable irrigation lines, 5 clips and 15 perforators.

**INCLUSIVE >**  
**1 Kit Essential**



## Accessoires

Promotion valable only for  
ImplantCenter2



**Newtron®  
LED handpiece**  
Ref. F12609  
For Piezotome® 2  
and ImplantCenter™ 2



**Newtron®  
LED cord**  
Ref. F12830  
For Piezotome® 2  
and ImplantCenter™ 2

~~€1.811~~  
**€1.090**



**Universal irrigation lines Sterile x10**  
Ref. F59905 - Single-use with integrated cassette and flow control system.

~~€174~~  
**€135**  
as from  
4 pieces



**Universal II Kit\***  
Ref. F87539  
Tips BS6, LC1, LC2,  
CE1, CE3 et PZ1.  
~~€1.356,85~~ **€1.090**



**Essential II Kit\***  
Ref. F87528  
Tips BS1S II, BS4 II,  
SL1 II, SL2 II, SL3 II,  
LC2 II.  
~~€1.356,85~~ **€1.090**



**Crest Splitting II Kit\***  
Ref. F87567  
Tips CS1 II, CS2 II, CS3 II,  
CS4 II, CS5 II, CS6 II.  
~~€878,39~~ **€819**



**Crown Extension II Kit\***  
Ref. F87554  
Tips BS6 II, CE1 II, CE2 II,  
CE3 II.  
~~€683,08~~ **€619**

# servotome

Delivered with a single footswitch,  
a box with 10 electrodes,  
a conductive bracelet.  
Ref. F57210

~~€1.976,45~~


**€1.550**





## Electrodes

Abscess incision, exposure of impacted tooth, cutting of peri-coronary sac, hemostasis, exposure of crown lengthening, gingivoplasty, gingivectomy, frenectomy, tissue ablation before impression, stripping, etc.


### Fulguration & coagulation




### Coagulating incision



### Incision



### Excision





# air-n-go easy

## All-in-one polisher for easy procedures

- > **AIR-N-GO® easy** with KaVo\* adapter Ref. F10100
- > **AIR-N-GO® easy** with Sirona\* adapter Ref. F10101
- > **AIR-N-GO® easy** with W&H\* adapter Ref. F10102
- > **AIR-N-GO® easy** with Bien-Air\* adapter Ref. F10103
- > **AIR-N-GO® easy** with Midwest\* adapter Ref. F10105



with a convertible AIR-N-GO® easy handpiece, a Supra nozzle, 1 CLASSIC powder bottle neutral (250 g), 1 AIR-N-GO® easy maintenance kit (replacement o'rings, syringe of silicone grease with brushes, cannula, cleaning drill).

\*KaVo, Sirona, W&H, Bien-Air and Midwest are registered trademarks belonging respectively to the company of the same name.

# air-n-go perio

air-n-go easy + air-n-go perio > ~~€2.070,30~~  
**€1.519**



### Kit PERIO easy

with a AIR-N-GO® easy nozzle and a bottle of AIR-N-GO® PERIO glycine powder 160g of AIR-N-GO® PERIO glycine powder.  
Ref. F10120

Dental air-polishing handpiece - Class IIa medical device (GMED) - CE 0459.

Manufacturer: SATELEC® ACTEON® Group - France.

## AIR-N-GO® "CLASSIC" powder

4 bottles 250 g

- > Peppermint Ref. F10250 ~~€102,70~~
- > Lemon Ref. F10251 **€99,71**
- > Raspberry Ref. F10252 each
- > Cola Ref. F10253 **€74,90**
- > Neutral Ref. F10254 as from
- > Assortment Ref. F10255 4 pieces



## air-n-go perio

### AIR-N-GO® "PERIO" powder

Ref. 10270  
3 bottles 160 g



~~€126,97~~  
**€112**

# miniled standard

~~€1.095,56~~

**€819**

Ref. F02507

## ► Fast

Only 6 to 12 sec to polymerize 2mm composite.

## ► Power

1250mW/cm<sup>2</sup>  
(Ø7,5mm light guide).

2000mW/cm<sup>2</sup>  
(Ø5,5mm light guide - optional).

# Highly Adaptable


The alliance of power & safety

## ► Safety

Base-charger includes a radiometer for light power intensity check-up.

## ► Optimal result


3 modes to fit all clinical needs.



### FAST MODE

#### Fast and easy


Full power during 10 sec  
(bip at 5 sec) on regular depth.



### PULSE MODE

#### Ally rapidity & safety

10 flashes of 1 sec  
(bip after 5 flashes)  
with times off to  
preserve dental pulp.



### RAMPING MODE

#### No stress cracking

Progressive power  
increase during 20 sec  
(bip every 5 sec) for  
a complete and safe  
curing.



# Simply Powerful

Outperforms in  
power, efficiency and speed

~~€1.407,36~~

**€1.090**

Ref. F05214



# miniled supercharged

## ► Visibility

Luminous screen with  
customizable modes and  
times, count curing time  
and load indicator.




## ► Power

2000mW/cm<sup>2</sup>  
(Ø7,5mm light guide).

3000mW/cm<sup>2</sup>  
(Ø5,5mm light guide - optional).

## ► Optimal result


3 modes to fit all clinical needs.



### FAST MODE

#### The choice of rapidity


Full power during 3, 4, 5 or 10  
sec, according to the composite  
thickness.



### PULSE MODE

#### Limits temperature rise

5 or 10 flashes of 1 sec for a  
polymerization in security.



### RAMPING MODE

#### Avoids composites retraction

Progressive power increase during  
9 or 15 sec to ensure perfect sealing  
of restoration.



# miniled ortho 2

~~€1.457,10~~

**€1.090**

Ref. F05220

## ► Customizable signal

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

## ► Power

3000mW/cm<sup>2</sup>  
(Ø5,5mm light guide).

## ► Long lasting battery

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

## ► Ortho friendly

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

## ► Efficacy & safety

Cure any type of material used in orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.

## Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from  
few brackets to a whole quadrant



# Mouth mirrors

## PURE REFLECT®

- > Very bright image
- > Resistant

For all purposes

### RHODIUM

- > A greyish reflection



### TITANIUM

- > A blueish reflection

## FRONT SURFACE Cross-sectional view






Image reflected directly on the surface of the glass, distortion-free  
**NO DOUBLE IMAGE**

## BACK SURFACE Cross-sectional view of an aluminium mirror



Reflecting surface under the glass  
**DOUBLE IMAGE**




## PURE REFLECT®

| N° | Ref.     | N. |
|----|----------|----|
| 4  | MP3622PH | 12 |
| 5  | MP3624PH | 12 |

~~€52,78~~

**€36**



## ALUMINIUM

| N° | Ref.     | N. |
|----|----------|----|
| 4  | MP3122PH | 12 |
| 5  | MP3124PH | 12 |

~~€18,72~~

**€13**



# MTA<sup>\*</sup>caps

## Endodontic repair cement

1 aluminium bag with 2 pre-proportioned 0,3g-capsules

Ref. 251170

~~€98,63~~ **€79**

# salvizol E.D.T.A

## Root canal irrigant

250ml bottle

Ref. 266115

~~€118,97~~ **€79**



# sealite regular

266264

## Powder and liquid for permanent sealing cement for root canals

25g of Sealite™ Regular Powder (1),  
10ml of Sealite™ Liquide (1), 1 spoon (2),  
1 dropper (2)

Ref. 266264

~~€123,62~~ **€79**

# desocclusol orange



## Solvent for eugenate root canal filling cements

Ref. 226022

~~€37,17~~ **€35**

|                                                                                      |                                                        |               |                   |            |
|--------------------------------------------------------------------------------------|--------------------------------------------------------|---------------|-------------------|------------|
|  | Pure Reflect® Micromirror Ø3mm - 14,8cm - round        | Ref. MP1303PH | <del>€75,71</del> | <b>€69</b> |
|  | Pure Reflect® Micromirror Ø5mm - 14,8cm - round        | Ref. MP1305PH | <del>€75,71</del> | <b>€69</b> |
|  | Pure Reflect® Micromirror 7x2mm - 14,8cm - rectangular | Ref. MP1365PH | <del>€87,17</del> | <b>€79</b> |
|  | Pure Reflect® Micromirror 9x3mm - 14,8cm - rectangular | Ref. MP1310PH | <del>€87,17</del> | <b>€79</b> |

# Riskontrol® contributes to clinical protection



## RISKONTROL® is also a GREEN solution

- > Disposable tip
- > Economic: simplify process to save time
- > Recyclable plastic: environmentally friendly

## Sterilizable tips are not sufficient anymore



## RISKONTROL® can be used for all procedures in any workflow




Class I and Class II are medical devices. Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of this device is in the Product Feature Summary available on the laboratory's website.

- > Cost-effective
- > Flexible for a better use: respect of patient's anatomy
- > Safety lock system: no risk of ejecting the tips
- > Astringent effect – no bleeding



## Use Expasyl™ Exact to improve patient comfort:

- > 37 less pressure than traditional cords
- > Atraumatic sulcular opening
- > Preserved epithelial attachment
- > No damage of the sulcus
- > Astringent effect - no bleeding



Expasyl™ generates pressure that is 1.7 to 9.2 times greater than other retraction paste available on the market.



# total protec by riskontrol

~~€32,12~~  
**€28**



**Box of 100 non-sterile tips with soft protection.**

- > White Ref. 201751
- > Blue Ref. 201750

# classic by riskontrol



~~€55,26~~  
**€31**

**Box of non-sterile tips.**

- > White Ref. 201664
- > Blue Ref. 201665
- > Yellow Ref. 201666
- > Green Ref. 201667
- > Purple Ref. 201668

# expasyl exact

~~€140,81~~  
**€95**

261010



**Paste for temporary sulcular opening  
20 capsules (0,3g) Expasyl™ Exact.**

# expasyl



~~€238,04~~  
**€144**

261030

**Paste for temporary sulcular opening  
20 capsules (1g) Expasyl™.**

## The X-Mind<sup>®</sup> generators can accomodate any operatory configuration.

3 arms lengths are available : 41 cm (16.2"), 82.5 cm (32.5") and 110 cm (43.5") and can be:

- > Top wall mounted
- > Bottom wall mounted

Dental unit, ceiling and mobile version are also available.

With its adaptable wall plate, it is easy to replace an older generator with the X-Mind<sup>®</sup> without having to do additional drilling on the wall.



~~€4.611,95~~

**€2.999**



| Wall version    | Top mounting    | Top mounting square cone | Bottom mounting | Bottom mounting square cone |
|-----------------|-----------------|--------------------------|-----------------|-----------------------------|
| 0,40 m extensor | Ref. 5.92.23137 | Ref. 5.92.23167          | Ref. 5.92.23130 | Ref. 5.92.23160             |
| 0,80 m extensor | Ref. 5.92.23138 | Ref. 5.92.23168          | Ref. 5.92.23131 | Ref. 5.92.23161             |
| 1,10 m extensor | Ref. 5.92.23139 | Ref. 5.92.23169          | Ref. 5.92.23132 | Ref. 5.92.23162             |


## Choose from a variety of cones to better fit to your clinical needs



### Long cone

(31 cm - 12")

Recommended for the parallel technique



### Short cone

(20 cm - 8")

Recommended for bisector technique



### Rectangular Long cone

(31 cm - 12")

Reduces the X-rays applied to the patient's skin by 50%.

# x-mind unity

## A sharp and contrasted image

- > A contrasted and defined image
- > Optimized protection for the patient
- > The generator and sensor association

Fixing from above, arm 0.80 m, short cone.  
Other configurations: see public price 2023.



## Intraoral X-ray

- > Very high quality
- > Minimum exposure time: 0,02s
- > Precise and flexible positioning
- > Programmable timer
- > Fast selection from traditional to digital mode
- > Top or bottom mounting

~~€5.626,40~~  
**€3.450**



X-MIND® DC et X-MIND® unity: Medical devices class IIb  
CE 0051 (IMQ). Manufacturer: de Götzen Groupe ACTEON® (Italy)

# sopix 2



## A remarkable image quality, even for endodontics and the more complexes

SOPIX<sup>2</sup> provides a high resolution image by exploiting the latest technologies available on the market.

### SOPIX<sup>2</sup> USB Sensor size 1

high definition with a sensor holder,  
10 sensor sheaths  
Ref S\_802\_0008

~~€6.751~~ **€3.990**

Also available in size 2  
Ref S\_802\_0009

~~€7.200~~ **€5.490**



# x-mind prime 2




~~€26.577,21~~  
**€15.500**

## X-MIND® Prime2 2D wall mounted

The ACTEON®'s smallest extraoral 2D imaging solution ever for the **maximum accessibility in any context**.

Ref. W1200001




~~€75.757,45~~  
**€43.690**

## X-MIND® Prime2 3D wall mounted

**Cutting-edge CBCT technology** enclosed in an essential and elegant style, fitting almost any space, even the narrowest ones.

Ref. W1200005


## Last detection technologies for high technical performance



**High resolution image,**  
up to 72.5µm of voxel size



**SMAR-T option,**  
the innovative Selective Metal  
Artifacts Reduction Tool



**Fast acquisition**



**FOV**

50x50, 85x93, 85x50, **120x100**

X-MIND® prime and X-MIND® prime 3D are manufactured by de Götzen (a company of the ACTEON® group) and distributed by ACTEON®. They are Class I medical devices with type B applied parts according to the IEC 60601-1 classification.





## X-MIND® Prime2 2D ceph

Cephalometric add-on for a full 2D system still encompassing the same philosophy of prime, **ideal for pediatric applications.**



Ref. W1200009



## X-MIND® Prime2 3D ceph

**The most comprehensive all-in-one ACTEON®** imaging solution, all you need to successfully handle any patient type in any clinical landscape... always at your fingertips.

Ref. W1200010



- > **High resolution image**, up to 72,5µm of voxel size
- > **SMAR-T**, innovative Selective Metal Artifacts Reduction Tool
- > **Field of view (FOV)**: 5x5 till 12x10

- > Compact design, **quick installation**
- > **Easy to use** software AIS 3D App
- > **Prime2 3D**, with pre-installed pc

# pspix2



—€8.154,19—  
**€ 6.490**



reddot award 2015  
winner



## The first personal imaging plate scanner

- > Compact design and easy to access: the smallest imaging plate scanner on the market for your personal convenience.
- > Images in seconds: you no longer need to move to another room to scan your images.
- > Featuring an exclusive “Click & Scan” concept, the new PSPIX® is totally intuitive: make your selection on the large touchscreen, drop in your imaging plate and let the PSPIX® do the rest...

Delivered with 4 ERLM (2x size 1, 2x size 2), 200 protective bags (100x size 1, 100x size 2), an imaging plate storage box, 1 stylus for touchscreen, 1 protection bag for PSPIX² system, 1 power supply 100-240 VAC, 1 Ethernet cable, Sopro® Imaging and his license. - Ref. S\_702\_0013




Class I Medical devices. Manufacturer: Sopro® - France.

## 2 FIBER PIXEL

## DIFFERENTIATION OF THE DENTAL TISSUE


FIBER2PIXEL® technology is based on the use of broad spectrum optical microfibres for the guided transmission of photon emissions in order to provide highly contrasted images.

Periodontics size 0    Periodontics size 1    Endodontics size 1    Orthodontics size 1




- > With FIBER2PIXEL® technology, PSPIX®² provides accurate images and striking contrast to ensure a reliable diagnosis.
- > Various sizes are available depending on patient morphology and clinical applications.


Cariology size 2



Implantology size 2



Bitewings size 3



## Protective bag for ERLM\*

For **PSPIX®<sup>2</sup>** and **PSPIX®**  
(previous generation)

~~€82,09~~  
**€56**



|        | <b>PSPIX®</b><br>(previous generation)                                            |         |
|--------|-----------------------------------------------------------------------------------|---------|
|        | Box of 250                                                                        |         |
| Size 0 |  | 700338L |
| Size 1 |  | 700339L |
| Size 2 |  | 700340L |
| Size 3 |  | 700341L |

## Capteurs ERLM\*

For **PSPIX®<sup>2</sup>**

~~€125,45~~  
**€89**



|                       | <b>PSPIX®<sup>2</sup></b>                                                           |        |
|-----------------------|-------------------------------------------------------------------------------------|--------|
|                       | Box of 2                                                                            |        |
| Size 0<br>(22 x 35mm) |  | 990215 |
| Size 1<br>(24 x 40mm) |  | 990216 |
| Size 2<br>(31 x 41mm) |  | 990217 |
| Size 3<br>(27 x 54mm) |  | 990218 |

## Protections for ERLM\*

~~€45,97~~  
**€25**

Protective cover for  
**PSPIX** and **PSPIX<sup>2</sup>**  
(previous generation)



| Compatible with             | <b>PSPIX<sup>2</sup></b> | <b>PSPIX</b> |
|-----------------------------|--------------------------|--------------|
| PSPIX®                      | ✓                        | ✓            |
| SOREDEX®<br>Digora          | ✓                        | ✓            |
| INSTRUMENTARIUM®<br>Express | ✓                        | ✓            |
| DÜRR DENTAL®<br>Vistascan®  | ✓                        | ✓            |
| KAVO®                       | ✓                        | ✓            |

|        | <b>PSPIX®</b><br>(previous generation)                                              |         |
|--------|-------------------------------------------------------------------------------------|---------|
|        | Box of 250                                                                          |         |
| Size 0 |  | 700342L |
| Size 1 |  | 700343L |
| Size 2 |  | 700344L |
| Size 3 |  | 700345L |

## Protective bag for ERLM\*

For **PSPIX®<sup>2</sup>**

~~€96,47~~  
**€63**



|        | <b>PSPIX®<sup>2</sup></b>                                                             |        |
|--------|---------------------------------------------------------------------------------------|--------|
|        | Box of 300                                                                            |        |
| Size 0 |  | 700546 |
| Size 1 |  | 700547 |
| Size 2 |  | 700548 |
| Size 3 |  | 700549 |

- > Protection against scratches, light and cross contamination
- > Ergonomic and easy to use
- > Economic solution

\* Luminescent Radio Memory Display for digital phosphor plate systems  
ERLM PSPIX<sup>2</sup> - Class IIa Medical devices (GMED) - CE0459.  
ERLM PSPIX - Class IIa Medical devices (Eurofins) - CE0537.  
Manufacturer: Sopro® (France)



# C50




~~€3.086~~

**€2.990**

(daylight)


Ref. S\_C50\_0002 DAY LIGHT

Daylight



Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.


Daylight<sup>+</sup>



**High contrasts** emphasizing anatomic fissures and details of the dentine structure / depth.

## 5 Visions

Macro



Intraoral



Tooth




Smile






## 4 Modes

### Daylight




Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

### Perio




Enhancement of **dental plaque and gingival inflammation** thanks to chromatic amplification.

### Daylight\*




**High contrasts** emphasizing anatomic fissures and details of the dentine structure/depth.

### Cario



Caries **detection by fluorescence technique** by LEDs.

### Portrait



- > Full HD: 1920x1080
- > Ergonomic design with small tip (easy to use with small mouth/children)
- > Autofocus
- > Diagnostic feature: dental plaque and caries detection (C50 Fluo)



~~€4.590~~

**€3.990**

(fluorescence technique)

Ref. S\_C50\_0001 FLUO

# Prophylaxis Workflow



**C50**

Ref. S\_C50\_0001 FLUO

Package deal:

~~€8.354,87~~

**€7.490**




**air-n-go  
easy**



**+ air-n-go  
perio**

Ref. F10120


**newtron  
p5XS**



Ref. F62100



**1**  
F00246




**1S**  
F00245



**H3**  
F00369



**H4L**  
F00114



**H4R**  
F00115