

VITALI V8 TOUCH



Dental units designed produced and tested in Italy



V8 TOUCH: THE TOP DENTAL UNIT OF THE VITALI PRODUCT RANGE



THE VITALI PRODUCT RANGE

An optimized and versatile product range - 3 models and multiple configurations to meet any customization need.



V8 TOUCH



T5 EVO PLUS



T5 MASTER

A more advanced choice - Each model offers the most advanced standards of design, functionality, operational ergonomics and comfort for the patient.

V8 TOUCH

ALL THE BEST
OF VITALI'S OFFER

+

THE PERFORMANCES
OF A TOUCH INTERFACE



TOUCH INTERFACE

THE BENEFITS OF INNOVATION

1. Simplicity of use: nowadays, almost everyone uses touch devices (smartphones, tablets, ticket reservations, ATMs).
2. Immediate understanding, thanks to the intuitiveness of the interface and of its functions.
3. Accuracy in parameter settings.
4. Swift management.



In summary: **WITH V8 TOUCH ALL OPERATIONS ARE MANAGED IN THE MOST SIMPLE AND EFFICIENT WAY**

2

TOP EFFICIENCY

V8 TOUCH IS THE DENTIST'S FIRST ASSISTANT

1. Up to 5 profiles of use can be stored for 5 different operators.
2. From each profile, you can manage customised settings, setting specific parameters of use for each instrument and for each type of treatment.
3. The settings can be recalled with a single command, which can also be a foot control.



In summary: **HOW TO HAVE 5 UNITS IN A SINGLE V8 TOUCH**

3

GREAT VERSATILITY

TO CHOOSE THE UNIT THAT BEST MATCHES
YOUR WORKING STYLE

1. **Continental** table with instrument arms.
2. **International** table with pending cords.
3. Version with a **Cart**.
4. **Bracket** or **column** cuspidor.
5. NEW nurse's console on pantographed arm or on arm with double articulation, with touch or 2-button keyboard.
6. The version with bracket cuspidor may be configured for left-handed users.



In summary: **THE SAME FUNCTIONALITY IN ALL THE VERSIONS**

4

NEW DESIGN

A STYLE THAT IS THE RESULT OF ITALIAN TOUCH

1. A unit with clean, contemporary, harmonious lines.
2. With rounded and welcoming shapes that ensure that the patient feels completely at ease.
3. Wide choice of colours for upholstery.



In summary: **A STYLE IN STEP WITH THE TIMES AND WITH TECHNOLOGY**

5

THE NEW VITALI CHAIR

MORE COMFORT FOR THE PATIENT

1. New chair that adapts to the patient's build.
2. Seat/backrest coordinated movement.
3. Anatomical headrest with pneumatic adjustment on 3 axes.
4. Wide choice of colours for upholstery.
5. On request: self-modelling Memory Foam padding and movable right armrest.



In summary: **EFFICIENCY IN THE NAME OF TOP COMFORT**




ERGONOMIC DESIGN

MORE FREEDOM OF MOVEMENT TO WORK BETTER

1. LIMITED SIZE AND LEANER CHAIR
2. WIDE RANGE OF MOVEMENT OF THE INSTRUMENT TABLE FOR EXCELLENT INSTRUMENT POSITIONING
3. 80 CM EXTENSION OF INSTRUMENT TUBES
4. 180° CUSPIDOR ROTATION
5. FOOT CONTROL WITHOUT CONNECTION CABLES

allow you to:

- find an ergonomic, more natural and less tiring, working position,
- optimize any type of intervention both when operating individually and when interacting with the nurse



In summary: **LESS FATIGUE AND MORE EFFICIENCY**



THE BEST PROTECTION TO PREVENT AND FIGHT ANY RISK

Mechanical safety devices

V8 TOUCH features anti-crushing safety devices for the chair backrest and seat, nurse's console and cuspidor.

Hygiene safety devices

1. Anti-reflux device
2. Removable and autoclavable details (instrument support in silicone, lamp handles, cuspidor accessories, tray cover)
3. Tray-holder in removable steel that can be sterilized
4. Removable chair paddings for easy and deep sanitation.
5. VDS Vitali Disinfection System (optional)




In summary: **SAFETY FOR ALL,
FOR PATIENTS AND THE DENTAL TEAM**



VITALI STRUCTURAL SOLIDITY

A HIDDEN VALUE THAT IS WORTHY IN TIME

1. The structural features of V8 TOUCH mirror the principles of Vitali's production philosophy: **structural solidity, resistance, reliability.**
2. The aluminium structure of the frame ensures resistance, lightness, perfect stability. Without any surprises in time.
3. The use of noble materials for the other details (including stainless steel, bronze, ceramics, special plastics) ensures resistance to any type of mechanical and chemical stress.



In summary: **STABLE AND ROBUST UNITS FOR INTENSIVE USE**

9

THE BEST LIGHT

V8 TOUCH USES ONLY HIGH-QUALITY OPERATING LAMPS

1. Choice between **Maia** and **Alya**.
2. LED lamps for maximum lighting in the operating field.
3. High-quality lamps for a uniform and well defined vision.
4. Perfect handling, great orientability.



In summary: **WORKING WITH HIGH-DEFINITION LAMPS**

10

ACCESSORIES AVAILABLE ON REQUEST TO CREATE YOUR OWN V8 TOUCH

1. **VDS certified system**


Automatic sanitizing of internal circuits that combines different disinfecting actions to ensure top safety (UV-C rays + Multiclean liquids).

2. **Multimedia suite**

Diagnostic data that are immediately available and that can be shared with the patient, ensuring an uninterrupted work flow.

3. **Memory Foam chair padding**

Thermosensitive material that self models itself to the patient for his comfort (no compression points).



In summary: **THE FULL OPTIONAL VERSION**


11

VITALI DISINFECTION SYSTEM THE MAXIMUM SAFETY OPTION

Each Vitali unit can have the VDS (Vitali Disinfection System) equipment: certified automatic and exclusive system for the effective disinfection of the inner circuits.

VDS can:

- ✓ treat the mains water with UV-C rays before being used by the dental unit (UV-Osmo-Eco);
- ✓ treat, with a wide-ranging disinfectant (Multiclean), the water circuit of operational instruments, reducing the risks of cross infection;
- ✓ use a sterile liquid as an alternative to the mains water to prevent the contamination of the water circuit;
- ✓ sterilize air intake;
- ✓ constantly monitor the functional parameters of the system to ensure their proper operation and maximum safety for the operator and the patient.



In summary: AN INCREASINGLY NECESSARY OPTION

WHY SHOULD YOU CHOOSE A VITALI UNIT?



V8 TOUCH

IT IS MADE IN ITALY

Italy is the world's third largest producer of technologies and furnishings for the dental sector, after the US and Germany.

About 60% of national production goes overseas.

In the dental sector, Italian technology enjoys international credibility.

MADE IN ITALY, MADE IN VITALI

- Vitali was founded more than 60 years ago and is a completely Italian company based near Bologna, in the centre of a district that is highly specialized in dental equipment and bio-medical devices.
- In practical terms, this means:
 - product reliability
 - ongoing services

VITALI INNOVATION

THE FACTS

VITALI has made important industrial investments

VITALI has invested in products

VITALI has strengthened its assistance and advice network in Italy and all over the world

THE EFFECTS

- **NEW PLANT**
- **FULL PRODUCTION CYCLE (design, production, control)**
- **RESPECT FOR THE ENVIRONMENT**

- **ENTIRELY NEW RANGE OF UNITS**

- **EFFICIENT SERVICE**
- **CUSTOMER SATISFACTION**

VITALI? NO PROBLEM!

How can you calculate the cost and value of a dental unit?
Only according to its purchase price? That's wrong.
You need to take into account also maintenance costs.

In the case of Vitali units, the low incidence of assistance interventions

has a double positive effect:

- **more profitability of the dental unit;**
- **lower total cost of ownership** (i.e. the cost that sums price + maintenance costs + any loss of income for inefficiency of the dental unit).

VITALI MADE FOR YOU

- Because it puts the professional at the heart of every project.
- Because it offers the right value for money ratio:
some solutions that can be used across the entire range of products, have allowed the company to optimize production costs, without sacrificing quality.



VITALI SOCIAL RESPONSIBILITY

- Vitali has been participating for several years in non profit projects designed to bring dental care to Eastern Europe countries.
- In particular, Vitali has endorsed the OVERLAND FOR SMILE onlus humanitarian project that aims at providing concrete assistance for the dental health of children in orphanages in Eastern Europe countries.
- 3 new units donated by Vitali to OVERLAND FOR SMILE, now equip a mobile dental studio installed on an Iveco truck.



WHAT THE DENTISTS SAY OF VITALI

"Top Functionality: and it is *indestructible* when it comes to reliability".

"Never had any problems".

"Aesthetically it is always as new: the plastics are perfect,
no colour change, no creaking".

"Had a small problem after many years: it was solved immediately".

"I have been working with them for years."

"It still has modern lines, even after several years:
you do *not get tired* of its design".

"You do not need assistance, but if you call them, they *always solve* any problems".



designed for reliability